

For Research Use Only

CXCR1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83801-2-PBS

Basic Information

Catalog Number:

83801-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_000634

GeneID (NCBI):

3577

UNIPROT ID:

P25024

Full Name:

interleukin 8 receptor, alpha

Calculated MW:

40 kDa

Purification Method:

Protein A purification

CloneNo.:

240819H5

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00731-1, CXCR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83801-2-PBS. Detection antibody: 83801-1-PBS. Standard: Ag28962. Range: 1.25-80 ng/mL



Cytometric bead array standard curve of MP00731-2, CXCR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83801-2-PBS. Detection antibody: 83801-3-PBS. Standard: Ag28962. Range: 5-80 ng/mL