

For Research Use Only

WDR21C Recombinant antibody, PBS Only (Detector)



Catalog Number: 83802-2-PBS

Basic Information

Catalog Number: 83802-2-PBS	GenBank Accession Number: BC111778	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 138009	CloneNo.: 240402F3
Source: Rabbit	UNIPROT ID: Q8NA75	
Isotype: IgG	Full Name: WD repeat domain 21C	
Immunogen Catalog Number: AG30748	Calculated MW: 395 aa, 44 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00738-1, WDR21C Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83802-3-PBS. Detection antibody: 83802-2-PBS. Standard: Ag30748. Range: 0.625-80 ng/mL.