

For Research Use Only

MITF Recombinant antibody, PBS Only (Capture)

Catalog Number: 83803-4-PBS



Basic Information

Catalog Number: 83803-4-PBS	GenBank Accession Number: BC012503	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4286	CloneNo.: 240815G9
Source: Rabbit	UNIPROT ID: O75030	
Isotype: IgG	Full Name: microphthalmia-associated transcription factor	
Immunogen Catalog Number: AG3679	Calculated MW: 91 aa, 10 kDa, 59 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00746-2, MITF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83803-4-PBS. Detection antibody: 83803-1-PBS. Standard: Ag3679. Range: 0.781-100 ng/mL.