

For Research Use Only

TBX21, T-bet Recombinant antibody, PBS Only (Capture)



Catalog Number: 83804-1-PBS

Basic Information

Catalog Number: 83804-1-PBS	GenBank Accession Number: BC039739	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 30009	CloneNo.: 240849A11
Source: Rabbit	ENSEMBL Gene ID: ENSG00000073861	
Isotype: IgG	UNIPROT ID: Q9UL17	
Immunogen Catalog Number: AG4618	Full Name: T-box 21	
	Calculated MW: 535 aa, 58 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00748-2, TBX21,T-bet Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83804-1-PBS. Detection antibody: 83804-3-PBS. Standard: Ag4618. Range: 0.313-40 ng/mL.