

For Research Use Only

# TBX21, T-bet Recombinant antibody, PBS Only (Detector)



Catalog Number: 83804-3-PBS

## Basic Information

Catalog Number: 83804-3-PBS	GenBank Accession Number: BC039739	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 30009	CloneNo.: 240849B7
Source: Rabbit	ENSEMBL Gene ID: ENSG00000073861	
Isotype: IgG	UNIPROT ID: Q9UL17	
Immunogen Catalog Number: AG4618	Full Name: T-box 21	
	Calculated MW: 535 aa, 58 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00748-2, TBX21,T-bet Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83804-1-PBS. Detection antibody: 83804-3-PBS. Standard: Ag4618. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00748-3, TBX21,T-bet Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83804-2-PBS. Detection antibody: 83804-3-PBS. Standard: Ag4618. Range: 0.313-40 ng/mL