

For Research Use Only

# CSNK2A1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83805-1-PBS

## Basic Information

|                                     |  |  |
|-------------------------------------|--|--|
| Catalog Number:<br>83805-1-PBS      | GenBank Accession Number:<br>BC011668              | Purification Method:<br>Protein A purification |
| Size:<br>1 mg/ml                    | GeneID (NCBI):<br>1457                             | CloneNo.:<br>240862A4                          |
| Source:<br>Rabbit                   | UNIPROT ID:<br>P68400                              |  |
| Isotype:<br>IgG                     | Full Name:<br>casein kinase 2, alpha 1 polypeptide |  |
| Immunogen Catalog Number:<br>AG1452 | Calculated MW:<br>45 kDa                           |  |

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00749-2, CSNK2A1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83805-2-PBS. Detection antibody: 83805-1-PBS. Standard: Ag1452. Range: 3.13-100 ng/mL.