

For Research Use Only

CSNK2A1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83805-3-PBS

Basic Information

Catalog Number: 83805-3-PBS	GenBank Accession Number: BC011668	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1457	CloneNo.: 240862A9
Source: Rabbit	UNIPROT ID: P68400	
Isotype: IgG	Full Name: casein kinase 2, alpha 1 polypeptide	
Immunogen Catalog Number: AG1452	Calculated MW: 45 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00749-1, CSNK2A1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83805-3-PBS. Detection antibody: 83805-2-PBS. Standard: Ag1452. Range: 0.78-100 ng/mL.