

For Research Use Only

GATA3 Recombinant antibody, PBS Only



Catalog Number: 83808-1-PBS

Basic Information

Catalog Number: 83808-1-PBS	GenBank Accession Number: BC003070	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2625	CloneNo.: 240684A2
Source: Rabbit	UNIPROT ID: P23771	
Isotype: IgG	Full Name: GATA binding protein 3	
Immunogen Catalog Number: AG0174	Calculated MW: 443 aa, 48 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
human, mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83808-1-RR (GATA3 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83808-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83808-1-RR against Human GATA3 were performed. The affinity constant is 0.686 nM.