

For Research Use Only

GATA3 Recombinant antibody, PBS Only



Catalog Number: 83808-2-PBS

Basic Information

Catalog Number: 83808-2-PBS	GenBank Accession Number: BC003070	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2625	CloneNo.: 24068488
Source: Rabbit	UNIPROT ID: P23771	
Isotype: IgG	Full Name: GATA binding protein 3	
Immunogen Catalog Number: AG0174	Calculated MW: 443 aa, 48 kDa	

Applications

Tested Applications:
ELISA, IF/ICC
Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed MCF-7 cells using GATA3 antibody (83808-2-RR, Clone: 240684B8) at dilution of 1:250 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 83808-2-PBS in a different storage buffer formulation.