

For Research Use Only

# Phospho-Histone H2B (Ser14) Recombinant antibody, PBS Only



Catalog Number: 83811-1-PBS

## Basic Information

Catalog Number: 83811-1-PBS	GenBank Accession Number: BC005827	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 8349	CloneNo.: 240570A8
Source: Rabbit	UNIPROT ID: Q16778	
Isotype: IgG	Full Name: histone cluster 2, H2be	
	Calculated MW: 14 kDa	

## Applications

Tested Applications:  
WB, ELISA  
Species Specificity:  
human, mouse

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Non-treated NIH/3T3 cells and PMA treated NIH/3T3 cells were subjected to SDS PAGE followed by western blot with 83811-1-RR (Phospho-Histone H2B (Ser14)) at dilution of 1:1000 incubated at room temperature for 1.5 hours. The membrane was stripped and re-blotted with Alpha Tubulin (66031-1-Ig) antibody as a loading control. This data was developed using the same antibody clone with 83811-1-PBS in a different storage buffer formulation.