

For Research Use Only

# ATG8L Recombinant antibody, PBS Only



Catalog Number: 83817-1-PBS

## Basic Information

<b>Catalog Number:</b> 83817-1-PBS	<b>GenBank Accession Number:</b> BC009309	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 23710	<b>CloneNo.:</b> 240533F5
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q9HOR8	
<b>Isotype:</b> IgG	<b>Full Name:</b> GABA(A) receptor-associated protein like 1	
<b>Immunogen Catalog Number:</b> AG1473	<b>Calculated MW:</b> 14 kDa	

## Applications

**Tested Applications:**  
ELISA, WB

**Species Specificity:**  
rat, mouse, human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83817-1-RR (ATG8L antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83817-1-PBS in a different storage buffer formulation.