

For Research Use Only

# SLC18A3 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83819-2-PBS

## Basic Information

Catalog Number:

83819-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25153

GenBank Accession Number:

BC007765

GeneID (NCBI):

6572

UNIPROT ID:

Q16572

Full Name:

solute carrier family 18 (vesicular acetylcholine), member 3

Calculated MW:

532 aa, 57 kDa

Purification Method:

Protein A purification

CloneNo.:

240687D2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00780-2, SLC18A3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83819-2-PBS. Detection antibody: 83819-1-PBS. Standard: Ag25153. Range: 1.563-100 ng/mL