

For Research Use Only

# GADD45GIP1 Recombinant antibody, PBS Only (Detector)



Catalog Number:83823-1-PBS

## Basic Information

Catalog Number:

83823-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9299

GenBank Accession Number:

BC013027

GeneID (NCBI):

90480

UNIPROT ID:

Q8TAE8

Full Name:

growth arrest and DNA-damage-  
inducible, gamma interacting protein  
1

Calculated MW:

222 aa, 25 kDa

Purification Method:

Protein A purification

CloneNo.:

240947A2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00800-3, GADD45GIP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83823-2-PBS. Detection antibody: 83823-1-PBS. Standard: Ag9299. Range: 0.78-100 ng/mL.