

For Research Use Only

GADD45GIP1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83823-2-PBS

Basic Information

Catalog Number: 83823-2-PBS	GenBank Accession Number: BC013027	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 90480	CloneNo.: 240947A6
Source: Rabbit	UNIPROT ID: Q8TAE8	
Isotype: IgG	Full Name: growth arrest and DNA-damage- inducible, gamma interacting protein 1	
Immunogen Catalog Number: AG9299	Calculated MW: 222 aa, 25 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00800-2, GADD45GIP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83823-2-PBS. Detection antibody: 83823-4-PBS. Standard: Ag9299. Range: 0.78-100 ng/mL.



Cytometric bead array standard curve of MP00800-3, GADD45GIP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83823-2-PBS. Detection antibody: 83823-1-PBS. Standard: Ag9299. Range: 0.78-100 ng/mL.