

For Research Use Only

CLCA1 Recombinant antibody, PBS Only (Capture)

Catalog Number:83829-4-PBS



Basic Information

Catalog Number: 83829-4-PBS	GenBank Accession Number: BC156805	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1179	CloneNo.: 240941F4
Source: Rabbit	UNIPROT ID: A8K7I4	
Isotype: IgG	Full Name: chloride channel accessory 1	
Immunogen Catalog Number: AG18176	Calculated MW: 914 aa, 100 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00794-1, CLCA1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83829-4-PBS. Detection antibody: 83829-1-PBS. Standard: Ag18176. Range: 0.625-80 ng/mL



Bilayer interferometry (BLI) kinetic assay of 83829-4-PBS against Human CLCA1 was performed. The affinity constant is below 1 pM.