

For Research Use Only

PROSC Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83831-1-PBS

Basic Information

Catalog Number: 83831-1-PBS	GenBank Accession Number: BC012334	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 11212	CloneNo.: 240998A6
Source: Rabbit	UNIPROT ID: O94903	
Isotype: IgG	Full Name: proline synthetase co-transcribed homolog (bacterial)	
Immunogen Catalog Number: AG18369	Calculated MW: 275 aa, 30 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00796-2, PROSC Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83831-1-PBS. Detection antibody: 83831-3-PBS. Standard: Ag18369. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP00796-3, PROSC Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83831-4-PBS. Detection antibody: 83831-1-PBS. Standard: Ag18369. Range: 1.56-100 ng/mL.