

For Research Use Only

# NHE3 Recombinant antibody, PBS Only (Capture)

Catalog Number:83832-5-PBS



## Basic Information

<b>Catalog Number:</b> 83832-5-PBS	<b>GenBank Accession Number:</b> BC101669	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 6550	<b>CloneNo.:</b> 240999F2
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P48764	
<b>Isotype:</b> IgG	<b>Full Name:</b> solute carrier family 9 (sodium/hydrogen exchanger), member 3	
<b>Immunogen Catalog Number:</b> AG25637	<b>Calculated MW:</b> 834 aa, 93 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00818-2, SLC9A3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83832-5-PBS. Detection antibody: 83832-3-PBS. Standard: Ag25637. Range: 0.625-40 ng/mL



Biolayer interferometry (BLI) kinetic assay of 83832-5-PBS against Human NHE3 was performed. The affinity constant is 1.00 nM.