

For Research Use Only

PTPRS Recombinant antibody, PBS Only (Capture)

Catalog Number: 83833-2-PBS



Basic Information

Catalog Number: 83833-2-PBS	GenBank Accession Number: NM_002850	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 5802	CloneNo.: 241005C5
Source: Rabbit	UNIPROT ID: Q13332	
Isotype: IgG	Full Name: protein tyrosine phosphatase, receptor type, S	
	Calculated MW: 217 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00820-3, rPTP α Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83833-2-PBS. Detection antibody: 83833-3-PBS. Standard: Eg0576. Range: 0.313-40 ng/mL