

For Research Use Only

DLL3 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83834-1-PBS

Basic Information

| | | |
|--|--|---|
| Catalog Number: 83834-1-PBS | GenBank Accession Number: BC000218 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 10683 | CloneNo.: 240956D3 |
| Source: Rabbit | UNIPROT ID: Q9NYJ7 | |
| Isotype: IgG | Full Name: delta-like 3 (Drosophila) | |
| Immunogen Catalog Number: AG0579 | Calculated MW: 65 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00798-1, DLL3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83834-2-PBS. Detection antibody: 83834-1-PBS. Standard: Ag0579. Range: 0.313-40 ng/mL