

For Research Use Only

DLL3 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83834-2-PBS

Basic Information

Catalog Number: 83834-2-PBS	GenBank Accession Number: BC000218	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10683	CloneNo.: 240956D5
Source: Rabbit	UNIPROT ID: Q9NYJ7	
Isotype: IgG	Full Name: delta-like 3 (Drosophila)	
Immunogen Catalog Number: AG0579	Calculated MW: 65 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00798-1, DLL3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83834-2-PBS. Detection antibody: 83834-1-PBS. Standard: Ag0579. Range: 0.313-40 ng/mL