

For Research Use Only

Beta Arrestin 1 Recombinant antibody, PBS Only (Capture)

Catalog Number:83837-3-PBS



Basic Information

Catalog Number: 83837-3-PBS	GenBank Accession Number: BC003636	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 408	CloneNo.: 240974F2
Source: Rabbit	UNIPROT ID: P49407	
Isotype: IgG	Full Name: arrestin, beta 1	
Immunogen Catalog Number: AG7608	Calculated MW: 47 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00811-1, Arrestin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83837-3-PBS. Detection antibody: 83837-4-PBS. Standard: Ag7608. Range: 0.625-80 ng/mL.