

For Research Use Only

TFPI Recombinant antibody, PBS Only (Detector)

Catalog Number:83838-2-PBS



Basic Information

Catalog Number:

83838-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_006287.6

GeneID (NCBI):

7035

UNIPROT ID:

P10646-1

Full Name:

tissue factor pathway inhibitor
(lipoprotein-associated coagulation inhibitor)

Calculated MW:

35 kDa

Purification Method:

Protein A purification

CloneNo.:

240977D5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00812-2, TFPI Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83838-1-PBS. Detection antibody: 83838-2-PBS. Standard: Eg0855. Range: 0.781-100 ng/mL.