

For Research Use Only

PRKCA Recombinant antibody, PBS Only (Capture)



Catalog Number: 83840-2-PBS

Basic Information

Catalog Number: 83840-2-PBS	GenBank Accession Number: BC109274	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 5578	CloneNo.: 240926B6
Source: Rabbit	UNIPROT ID: P17252	
Isotype: IgG	Full Name: protein kinase C, alpha	
Immunogen Catalog Number: AG1025		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00814-2, PRKCA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83840-2-PBS. Detection antibody: 83840-3-PBS. Standard: Ag1025. Range: 0.625-80 ng/mL