

For Research Use Only

SMAD2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83841-1-PBS

Basic Information

Catalog Number: 83841-1-PBS	GenBank Accession Number: BC014840	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4087	CloneNo.: 240950D4
Source: Rabbit	UNIPROT ID: Q15796	
Isotype: IgG	Full Name: SMAD family member 2	
Immunogen Catalog Number: AG3237	Calculated MW: 467 aa, 52 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00791-2, SMAD2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83841-3-PBS. Detection antibody: 83841-1-PBS. Standard: Ag3237. Range: 0.469-60 ng/mL



Cytometric bead array standard curve of MP00791-3, SMAD2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83841-2-PBS. Detection antibody: 83841-1-PBS. Standard: Ag3237. Range: 0.938-60 ng/mL