

For Research Use Only

SLC25A20 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83842-1-PBS

Basic Information

Catalog Number: 83842-1-PBS	GenBank Accession Number: BC001689	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 788	CloneNo.: 240976C4
Source: Rabbit	UNIPROT ID: O43772	
Isotype: IgG	Full Name: solute carrier family 25 (carnitine/acylcarnitine translocase), member 20	
Immunogen Catalog Number: AG6896	Calculated MW: 33 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00793-1, SLC25A20 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83842-2-PBS. Detection antibody: 83842-1-PBS. Standard: Ag6896. Range: 0.781-100 ng/mL.