

For Research Use Only

TRMT6 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83843-1-PBS

Basic Information

Catalog Number: 83843-1-PBS	GenBank Accession Number: BC001262	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 51605	CloneNo.: 240911A1
Source: Rabbit	UNIPROT ID: Q9UJA5	
Isotype: IgG	Full Name: tRNA methyltransferase 6 homolog (S. cerevisiae)	
Immunogen Catalog Number: AG10176	Calculated MW: 56 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00816-2, TRMT6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83843-1-PBS. Detection antibody: 83843-4-PBS. Standard: Ag10176. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00816-3, TRMT6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83843-1-PBS. Detection antibody: 83843-3-PBS. Standard: Ag10176. Range: 0.625-80 ng/mL