

For Research Use Only

SLC31A1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83844-4-PBS

Basic Information

Catalog Number: 83844-4-PBS	GenBank Accession Number: BC013611	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1317	CloneNo.: 240868F10
Source: Rabbit	UNIPROT ID: O15431	
Isotype: IgG	Full Name: solute carrier family 31 (copper transporters), member 1	
Immunogen Catalog Number: AG26601	Calculated MW: 21 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00817-2, SLC31A1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83844-2-PBS. Detection antibody: 83844-4-PBS. Standard: Ag26601. Range: 0.313-20 ng/mL.