

For Research Use Only

RPS17 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83848-4-PBS

Basic Information

| | | |
|-------------------------------------|---------------------------------------|--|
| Catalog Number: 83848-4-PBS | GenBank Accession Number: BC009407 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 6218 | CloneNo.: 240851H3 |
| Source: Rabbit | UNIPROT ID: P08708 | |
| Isotype: IgG | Full Name: ribosomal protein S17 | |
| Immunogen Catalog Number: AG9334 | Calculated MW: 135 aa, 16 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00785-2, RPS17 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83848-4-PBS. Detection antibody: 83848-3-PBS. Standard: Ag9334. Range: 1.25-40 ng/mL.