

For Research Use Only

# FRMD3 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83850-3-PBS

## Basic Information

Catalog Number: 83850-3-PBS	GenBank Accession Number: BC023560	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 257019	CloneNo.: 240967G1
Source: Rabbit	UNIPROT ID: A2A2Y4	
Isotype: IgG	Full Name: FERM domain containing 3	
Immunogen Catalog Number: AG13235	Calculated MW: 69 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00787-1, FRMD3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83850-3-PBS. Detection antibody: 83850-2-PBS. Standard: Ag13235. Range: 0.78-100 ng/mL



Cytometric bead array standard curve of MP00787-2, FRMD3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83850-3-PBS. Detection antibody: 83850-1-PBS. Standard: Ag13235. Range: 0.78-100 ng/mL