

For Research Use Only

# PTPN12 Recombinant antibody, PBS Only (Detector)



Catalog Number:83851-1-PBS

## Basic Information

Catalog Number:

83851-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35382

GenBank Accession Number:

NM\_002835.3

GeneID (NCBI):

5782

UNIPROT ID:

Q05209

Full Name:

protein tyrosine phosphatase, non-receptor type 12

Calculated MW:

88kDa

Purification Method:

Protein A purification

CloneNo.:

240983A2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00788-3, PTPN12 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83851-4-PBS. Detection antibody: 83851-1-PBS. Standard: Ag35382. Range: 0.47-60 ng/mL.