

For Research Use Only

PTPN12 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83851-5-PBS

Basic Information

Catalog Number: 83851-5-PBS	GenBank Accession Number: NM_002835.3	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 5782	CloneNo.: 240983E1
Source: Rabbit	UNIPROT ID: Q05209	
Isotype: IgG	Full Name: protein tyrosine phosphatase, non-receptor type 12	
Immunogen Catalog Number: AG35382	Calculated MW: 88kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00788-2, PTPN12 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83851-2-PBS. Detection antibody: 83851-5-PBS. Standard: Ag35382. Range: 0.47-60 ng/mL.