

For Research Use Only

vwf Recombinant antibody, PBS Only (Capture)

Catalog Number: 83854-3-PBS



Basic Information

Catalog Number: 83854-3-PBS	GenBank Accession Number: GeneID (NCBI): 7450	Purification Method: Protein A purification
Size: 1 mg/ml	UNIPROT ID: P04275	CloneNo.: 240867H2
Source: Rabbit	Full Name: von Willebrand factor	
Isotype: IgG		
Immunogen Catalog Number: AG25578		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00810-1, VWF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83854-3-PBS. Detection antibody: 83854-1-PBS. Standard: Ag25578. Range: 0.313-40 ng/mL