

For Research Use Only

SOD1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83857-5-PBS

Basic Information

Catalog Number: 83857-5-PBS	GenBank Accession Number: NM_000454.5	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 6647	CloneNo.: 240888H6
Source: Rabbit	UNIPROT ID: P00441	
Isotype: IgG	Full Name: superoxide dismutase 1, soluble	
	Calculated MW: 16kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00799-3, SOD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83857-5-PBS. Detection antibody: 83857-2-PBS. Standard: Eg1043. Range: 0.313-20 ng/mL