

For Research Use Only

# RPE65 Recombinant antibody, PBS Only



Catalog Number:83861-1-PBS

## Basic Information

Catalog Number:

83861-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12356

GenBank Accession Number:

BC075036

GeneID (NCBI):

6121

UNIPROT ID:

Q16518

Full Name:

retinal pigment epithelium-specific protein 65kDa

Calculated MW:

61 kDa

Purification Method:

Protein A purification

CloneNo.:

240518F1

## Applications

Tested Applications:

ELISA, WB

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83861-1-RR (RPE65 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83861-1-PBS in a different storage buffer formulation.