

For Research Use Only

FAM127B Recombinant antibody, PBS Only



Catalog Number:83862-2-PBS

Basic Information

Catalog Number: 83862-2-PBS	GenBank Accession Number: BC000393	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 26071	CloneNo.: 240716B7
Source: Rabbit	UNIPROT ID: Q9BWD3	
Isotype: IgG	Full Name: family with sequence similarity 127, member B	
Immunogen Catalog Number: AG13930	Calculated MW: 113 aa, 13 kDa	

Applications

Tested Applications:
IF/ICC, ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed HepG2 cells using FAM127B antibody (83862-2-RR, Clone: 240716B7) at dilution of 1:250 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 83862-2-PBS in a different storage buffer formulation.