

For Research Use Only

# IGF1R Recombinant antibody, PBS Only (Capture)

Catalog Number: 83873-5-PBS



## Basic Information

Catalog Number:

83873-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_000875

GeneID (NCBI):

3480

UNIPROT ID:

P08069

Full Name:

IGF I receptor

Calculated MW:

155 kDa

Purification Method:

Protein A purification

CloneNo.:

240801D10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00869-1, IGF1R Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83873-5-PBS. Detection antibody: 83873-4-PBS. Standard: Ag24732 . Range: 0.156-10 ng/mL