

For Research Use Only

CD180 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83874-2-PBS

Basic Information

Catalog Number: 83874-2-PBS	GenBank Accession Number: BC109069	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4064	CloneNo.: 240910C12
Source: Rabbit	UNIPROT ID: Q99467	
Isotype: IgG	Full Name: CD180 molecule	
Immunogen Catalog Number: AG16444	Calculated MW: 661 aa, 74 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00825-1, CD180 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83874-1-PBS. Detection antibody: 83874-2-PBS. Standard: Ag16444. Range: 2.5-80 ng/mL.