

For Research Use Only

PCDHGB3 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83879-1-PBS

Basic Information

Catalog Number: 83879-1-PBS	GenBank Accession Number: BC150172	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 56102	CloneNo.: 241052A3
Source: Rabbit	UNIPROT ID: Q9Y5G1	
Isotype: IgG	Full Name: protocadherin gamma subfamily B, 3	
Immunogen Catalog Number: AG15569	Calculated MW: 929 aa, 101 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00862-1, PCDHGB3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83879-2-PBS. Detection antibody: 83879-1-PBS. Standard: Ag15569. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP00862-3, PCDHGB3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83879-3-PBS. Detection antibody: 83879-1-PBS. Standard: Ag15569. Range: 0.313-40 ng/mL.