

For Research Use Only

NCR3/NKp30 Recombinant antibody, PBS Only (Detector)

Catalog Number:83882-1-PBS



Basic Information

Catalog Number: 83882-1-PBS	GenBank Accession Number: BC052582	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 259197	CloneNo.: 240809A9
Source: Rabbit	UNIPROT ID: O14931	
Isotype: IgG	Full Name: natural cytotoxicity triggering receptor 3	
	Calculated MW: 201 aa, 22 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00833-3, NCR3/NKp30 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83882-3-PBS. Detection antibody: 83882-1-PBS. Standard: Eg31663. Range: 0.313-40 ng/mL.