

For Research Use Only

NCR3/NKp30 Recombinant antibody, PBS Only (Detector)

Catalog Number:83882-2-PBS



Basic Information

Catalog Number:

83882-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC052582

GeneID (NCBI):

259197

UNIPROT ID:

O14931

Full Name:

natural cytotoxicity triggering
receptor 3

Calculated MW:

201 aa, 22 kDa

Purification Method:

Protein A purification

CloneNo.:

240809C7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00833-1, NCR3/NKp30 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83882-4-PBS. Detection antibody: 83882-2-PBS. Standard: Eg31663. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP00833-2, NCR3/NKp30 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83882-3-PBS. Detection antibody: 83882-2-PBS. Standard: Eg31663. Range: 0.313-40 ng/mL.