

For Research Use Only

NCR3/NKp30 Recombinant antibody, PBS Only (Detector)

Catalog Number:83882-5-PBS



Basic Information

Catalog Number: 83882-5-PBS	GenBank Accession Number: BC052582	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 259197	CloneNo.: 240809B11
Source: Rabbit	UNIPROT ID: O14931	
Isotype: IgG	Full Name: natural cytotoxicity triggering receptor 3	
	Calculated MW: 201 aa, 22 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00833-4, Human NCR3/NKp30 Recombinant Matched Antibody Pair - PBS only. 83882-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg31663. 83882-5-PBS was HRP conjugated as the detection antibody. Range: 15.6-1000 pg/mL



Serum of eight individual healthy human donors was measured. The NCR3 concentration of detected samples was determined to be 143.6 pg/mL with a range of 77.4 - 279.5 pg/mL.