

For Research Use Only

Mouse CD200R1 Recombinant antibody, PBS Only (Capture)

Catalog Number:83889-1-PBS



Basic Information

Catalog Number:

83889-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_021325.3

GeneID (NCBI):

57781

UNIPROT ID:

Q9ES57

Full Name:

CD200 receptor 1

Calculated MW:

36 kDa

Purification Method:

Protein A purification

CloneNo.:

240979A10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00835-1, MOUSE CD200R1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83889-1-PBS. Detection antibody: 83889-2-PBS. Standard: Eg0988 Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00835-2, MOUSE CD200R1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83889-1-PBS. Detection antibody: 83889-3-PBS. Standard: Eg0988. Range: 0.313-40 ng/mL