

For Research Use Only

4EBP1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83890-4-PBS

Basic Information

Catalog Number: 83890-4-PBS	GenBank Accession Number: BC004459	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1978	CloneNo.: 241055H1
Source: Rabbit	UNIPROT ID: Q13541	
Isotype: IgG	Full Name: eukaryotic translation initiation factor 4E binding protein 1	
Immunogen Catalog Number: AG5056	Calculated MW: 118 aa, 12 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00559-4, 4EBP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83890-4-PBS. Detection antibody: 83890-1-PBS. Standard: Ag5056. Range: 0.625-80 ng/mL.