

For Research Use Only

# PHF19 Recombinant antibody, PBS Only (Detector)



Catalog Number:83892-2-PBS

## Basic Information

Catalog Number: 83892-2-PBS	GenBank Accession Number: BC022374	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 26147	CloneNo.: 241047B11
Source: Rabbit	UNIPROT ID: Q5T6S3	
Isotype: IgG	Full Name: PHD finger protein 19	
Immunogen Catalog Number: AG2511	Calculated MW: 580 aa, 63 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00868-3, PHF19 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83892-6-PBS. Detection antibody: 83892-2-PBS. Standard: Ag2511 Range: 0.313-40 ng/mL