

For Research Use Only

MPZ / PO Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number:83895-2-PBS



Basic Information

Catalog Number: 83895-2-PBS	GenBank Accession Number: BC006491	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4359	CloneNo.: 240948C3
Source: Rabbit	UNIPROT ID: P25189	
Isotype: IgG	Full Name: myelin protein zero	
Immunogen Catalog Number: AG0848	Calculated MW: 28 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00834-1, MPZ / PO Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83895-3-PBS. Detection antibody: 83895-2-PBS. Standard: Ag0848. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP00834-2, MPZ / PO Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83895-2-PBS. Detection antibody: 83895-1-PBS. Standard: Ag0848. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP00834-3, MPZ,PO Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83895-4-PBS. Detection antibody: 83895-2-PBS. Standard: Ag0848. Range: 1.563-100 ng/mL