

For Research Use Only

SIK1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83897-2-PBS

Basic Information

Catalog Number: 83897-2-PBS	GenBank Accession Number: BC038504	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 150094	CloneNo.: 240961C4
Source: Rabbit	UNIPROT ID: P57059	
Isotype: IgG	Full Name: salt-inducible kinase 1	
Immunogen Catalog Number: AG28481	Calculated MW: 783 aa, 85 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00838-2, SIK1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83897-4-PBS. Detection antibody: 83897-2-PBS. Standard: Ag28481. Range: 0.625-40 ng/mL