

For Research Use Only

FOXN2 Recombinant antibody, PBS Only (Detector)



Catalog Number:83899-4-PBS

Basic Information

Catalog Number: 83899-4-PBS	GenBank Accession Number: BC063305	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3344	CloneNo.: 240896H4
Source: Rabbit	UNIPROT ID: P32314	
Isotype: IgG	Full Name: forkhead box N2	
Immunogen Catalog Number: AG5126	Calculated MW: 47 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00842-1, FOXN2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83899-1-PBS. Detection antibody: 83899-4-PBS. Standard: Ag5126. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP00842-2, FOXN2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83899-3-PBS. Detection antibody: 83899-4-PBS. Standard: Ag5126. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP00842-3, FOXN2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83899-2-PBS. Detection antibody: 83899-4-PBS. Standard: Ag5126. Range: 1.563-100 ng/mL.