

For Research Use Only

KIAA1586 Recombinant antibody, PBS Only (Detector)



Catalog Number:83900-3-PBS

Basic Information

Catalog Number: 83900-3-PBS	GenBank Accession Number: BC041125	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 57691	CloneNo.: 240972D7
Source: Rabbit	UNIPROT ID: Q9HCI6	
Isotype: IgG	Full Name: KIAA1586	
Immunogen Catalog Number: AG36038	Calculated MW: 787 aa, 90 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00847-2, KIAA1586 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83900-6-PBS. Detection antibody: 83900-3-PBS. Standard: Ag36038. Range: 0.781-100 ng/mL