

For Research Use Only

VPRBP Recombinant antibody, PBS Only (Detector)



Catalog Number:83901-4-PBS

Basic Information

Catalog Number: 83901-4-PBS	GenBank Accession Number: BC022792	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 9730	CloneNo.: 241028G4
Source: Rabbit	UNIPROT ID: Q9Y4B6	
Isotype: IgG	Full Name: Vpr (HIV-1) binding protein	
Immunogen Catalog Number: AG2184	Calculated MW: 1506 aa, 169 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00852-1, VPRBP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83901-1-PBS. Detection antibody: 83901-4-PBS. Standard: Ag2184. Range: 0.625-40 ng/mL



Cytometric bead array standard curve of MP00852-2, VPRBP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83901-3-PBS. Detection antibody: 83901-4-PBS. Standard: Ag2184. Range: 0.313-40 ng/mL