

For Research Use Only

NDRG2 Recombinant antibody, PBS Only (Capture)



Catalog Number:83902-1-PBS

Basic Information

Catalog Number: 83902-1-PBS	GenBank Accession Number: BC010458	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 57447	CloneNo.: 241039A6
Source: Rabbit	UNIPROT ID: Q9UN36	
Isotype: IgG	Full Name: NDRG family member 2	
Immunogen Catalog Number: AG2629	Calculated MW: 371 aa, 39 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00828-2, NDRG2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83902-1-PBS. Detection antibody: 83902-3-PBS. Standard: Ag2629. Range: 0.313-20 ng/mL